



Basic Endoscopy Kit

One part delivers everything you need.
From setup through cleanup, we have it covered.

When there's less time for procedure preparation and cleanup, leave it to Medovations to package everything you need conveniently. Our Basic Endoscopy Kit contains some of the fundamental items needed for upper and lower GI procedures. You save purchasing time and labor setup by pulling one simple kit from inventory.

Open the bundle to reveal a generously-sized transport pad that holds lubricating jelly, suction tubing, cotton gauze pads, a compact wash basin and a water-activated sponge imbibed with enzymatic cleanser.

After the procedure, dispose of everything and use the cinch cords to turn the transport pad into a secure pouch to easily transfer the endoscope for reprocessing.



Add water to the basin to properly dilute enzymatic cleanser in the CleanFreak® Dry Imbibed Sponge.



Cinch the cords on the transport pad to create a bag to carry the endoscope for reprocessing.

Basic Endoscopy Kit

Part No.	Description	Qty/Case
1201-50	<ul style="list-style-type: none"> • 1 Transport Pad • 1 Water Basin • 1 Dry Imbibed Endoscope Sponge • 1 Tube E-Z Lubricating Jelly • 10 Cotton Gauze Pads • 1 Suction Tubing 	12
1201-60	<ul style="list-style-type: none"> • 1 Transport Pad • 1 Water Basin • 1 Dry Imbibed Endoscope Sponge • 1 Tube E-Z Lubricating Jelly • 10 Cotton Gauze Pads • 1 CleanFreak® Combination Brush 	20

Learn more from your personal product representative at **800-558-6408** or **Medovations.com**



Medovations, Inc. | 102 East Keefe Avenue | Milwaukee, Wisconsin 53212 USA
800.558.6408 or 414.265.7620 | Fax 414.265.7628
Orders: orders@medovations.com | General Inquiries: medo@medovations.com