



BullDog® Biopsy Valves and Accessories

Presenting the BullDog® Biopsy Valve.



Accommodate larger endoscopic instruments by opening the cap.

Cleaning an endoscope properly after a procedure is a critical part of the process. But when cleaning simple components becomes time consuming and ineffective, look to Medovations for economical, single-use solutions.

Reduce the risk of cross contamination with a single-use valve.

Reusable biopsy valves are difficult to clean effectively, and not worth the time or risk. Our single-use biopsy valves provide a no-worry solution that maintains secure access for endoscopic instruments, and reduces the risk of cross contamination.

Engineered for a perfect fit.

The biopsy valve is designed to seal around the instrument as it is passed down the accessory channel of the endoscope. To accommodate larger instruments, lift the cap to reveal a larger entry port. Color coding ensures the appropriate valve is matched to the endoscope brand.



Utilize the Irrigating Adapter to flush and clear the view.



The valve's slit opening seals to help prevent leakage and back flow.

BullDog® Endoscope Biopsy Valve

Part No.	Description	Qty/Case	Color
1203-00	Olympus/Fujinon Compatible	100	Blue
1203-10	Olympus/Fujinon Compatible, Sterile	100	Blue
1203-20	Pentax Compatible	100	Green
1203-50	Irrigating Adapter	50	Purple

Single-use. Made in the USA.

Learn more from your personal product representative at **800-558-6408** or **Medovations.com**



Medovations, Inc. | 102 East Keefe Avenue | Milwaukee, Wisconsin 53212 USA
 800.558.6408 or 414.265.7620 | Fax 414.265.7628
 Orders: orders@medovations.com | General Inquiries: medo@medovations.com