



# Capture™ Polyp Trap



# Presenting the Capture™ Polyp Trap.

The smart solution to capturing, collecting and transporting samples.

## Capture™ makes everything easier.

At Medovations, we're committed to making products that help simplify complicated procedures. By combining top-of-the-line elements with a user-friendly design, the new Capture™ Polyp Trap does just that.

## Simplify the Process.

The Capture™ Polyp Trap offers a range of features that makes the process of capturing, collecting and transporting specimens easier. Use one product to collect more polyps with less effort.

## Effortless Capture.

The transparent, fixed chambers with built-in strainers feature large, easy-to-read numbers and unobstructed visibility.

When suctioning, simply align the integral lid tab with the numbered chamber to capture the specimen.

## Convenient Transport.

Connect the PVC tube to the suction connector to seal the Capture™ Polyp Trap for safe transport and handling post procedure.

## Capture™ Polyp Trap

Part No.	Description	Qty/Case
1202-01	Polyp Trap—4 Chamber	25

Single-use, non-sterile. Made in the USA.



Four chambers with built-in strainers collect and separate specimens.



Easy-to-read numbers label each chamber.

Learn more from your personal product representative at **800-558-6408** or **Medovations.com**



Medovations, Inc. | 102 East Keefe Avenue | Milwaukee, Wisconsin 53212 USA  
 800.558.6408 or 414.265.7620 | Fax 414.265.7628  
 Orders: [orders@medovations.com](mailto:orders@medovations.com) | General Inquiries: [medo@medovations.com](mailto:medo@medovations.com)