



CleanFreak® Endoscope Transport Pad



Introducing the CleanFreak[®] Endoscope Transport Pad.

Simplify transport from procedure to reprocessing.

The CleanFreak Endoscope Transport Pad is a convenient post-procedure solution for handling endoscopes and other accessories.

Convenient Design

The CleanFreak Endoscope Transport Pad quickly unfolds to 40" x 29" with a flexible drape design suitable for use on large and small surfaces.

Safe, Smart Transport

After the procedure, cinching the strings on both ends of the pad will convert it into a leak-proof bag for transporting equipment for reprocessing.

The plastic liner features absorbent padding to contain fluids, and is sized generously to ensure the endoscope is not coiled too tightly during transport.



Edge-to-edge absorbent padding



Pull cords around dirty endoscope to transport

CleanFreak[®] Endoscope Transport Pad

| Part No. | Size | Qty/Case |
|----------|-----------|----------|
| 1201-90 | 40" x 29" | 25 |

Single-use, non-sterile.

Supports SGNA Guidelines* for Infection Control

- Transport the soiled endoscope to the reprocessing area in a closed container
- Containers, sinks, and basin should be large enough that the endoscope will not be damaged by being coiled too tightly.

*Standards of Infection Control in Reprocessing of Flexible Gastrointestinal Endoscopes, SGNA © 2012

Learn more from your personal product representative at **800-558-6408** or **Medovations.com**



Medovations, Inc. | 102 East Keefe Avenue | Milwaukee, Wisconsin 53212 USA
 800.558.6408 or 414.265.7620 | Fax 414.265.7628
 Orders: orders@medovations.com | General Inquiries: medo@medovations.com