


Flexible Yankauer Suction Device



- Flexible version of the traditional Yankauer suction device
- Suction tip can be formed to a precise suctioning angle during procedure
- Flexible suction end minimizes tissue trauma during suctioning
- Thermoformed suction end returns to original angle when not in use
- Available in 3 volume capacities; fine, regular, and high
- Suction Control Port (available on regular capacity) gives end-user complete control of vacuum pressure
- Basket Tip Guard (available on regular capacity) prevents tissue and bone chips from occluding the device
- Sterile

Suction Control Port and Basket Tip Guard options

Part Number	Suction Capacity	Suction Control Port	Basket Tip Guard	Qty per cs
1283-01	Fine	No	No	50
1284-01	Regular	No	No	50
1284-02	Regular	No	Yes	50
1284-03	Regular	Yes	No	50
1284-04	Regular	Yes	Yes	50
1285-01	High	No	No	50