

# Infiltration Tubing



## Infiltration Tubing

**Model SU1161—Dual Spike Tubing**

**Model SU1163—Single Spike Tubing**

Large bore “memory-free” tubing in full eight foot lengths with a luer lock fitting on the operating end. Available in dual and single spike configurations.

## Infiltration Tubing for Peristaltic Pumps

**Model SU1168**

**Model SU1170-LY—Lysonix® compatible**

Infiltration infusion pump tubing has a nine-foot length and a luer lock fitting. Peristaltic silicone tubing is pressure-locked to a six-foot section of tubing with a bag spike. The silicone tubing section is designed for increased flow and easy insertion on the pump head.

Model	Description	Count Per Case
SU1161	Infiltration Dual Spike Tubing, sterile	10
SU1163	Infiltration Single Spike Tubing, sterile	10
SU1168	Infiltration Tubing for Peristaltic Pump, sterile	10
SU1170-LY	Bifurcated Infiltration Tubing for LySonix® Irrigator, sterile	10



SU1161



SU1163



SU1168



SU1170-LY



**medovations®**

Endoscopy • Liposuction • Neurosurgery

