


5-in-1 Connectors

- Stepped
- Smooth



STEPPED



- Inner diameter expands proportionately with outer diameter
- Polypropylene material can be easily cut to match ID of connector to ID of tube
- Cut connector maximizes flow rate and prevents fluid buildup at connection point
- Reduces inventory of connectors by accommodating five different tube sizes
- 5 connection points on 1 connector
- Sterile or non-sterile

SMOOTH



- Reduces inventory of connectors by accommodating five different tube sizes
- 5 connection points on 1 connector
- Non-sterile

Part Number	Description	Material	Fits Tubes	Qty per cs
<i>Stepped:</i>				
1350-00	Non-Sterile	Polypropylene	5/16" - 9/16" (8mm-14mm)	100
1350-55	Sterile	Polypropylene	5/16" - 9/16" (8mm-14mm)	100
<i>Smooth:</i>				
1350-75	Non-Sterile	Polypropylene	3/16" - 7/16" (5mm-11mm)	100

102 E. KEEFE AVENUE • MILWAUKEE, WISCONSIN 53212 U.S.A.

Tel 414/265•7620 • Fax 414/265•7628

e-mail: medo@medovations.com

1•800•558•6408