




Thoracic Catheter

- Straight
- Right Angle



- Clear, kink-resistant, medical grade PVC material
- Angle-cut end for easy insertion
- QuickConnect™ funnel end can be cropped to form integral funnel connector
- QuickConnect™ funnel connector provides quick, easy connection to standard drainage units
- Smooth drainage eyelets and open distal tip enhance efficient drainage and maximize patient comfort
- Radiopaque stripe
- Depth markings in 2 cm increments
- Right Angle style available
- Adult, pediatric, and infant sizes
- Sterile

Part Number	French Size	O.D. (mm)	Drainage Eyelets	Length	Qty per cs	Part Number	French Size	O.D. (mm)	Drainage Eyelets	Length	Qty per cs
STRAIGHT:						RIGHT ANGLE:					
1240-12	12	4.0	4	20"(51cm)	10	1245-16	16	5.3	4	20"(51cm)	10
1240-16	16	5.3	4	20"(51cm)	10	1245-20	20	6.7	5	20"(51cm)	10
1240-20	20	6.7	5	20"(51cm)	10	1245-24	24	8.0	6	20"(51cm)	10
1240-24	24	8.0	6	20"(51cm)	10	1245-28	28	9.3	6	20"(51cm)	10
1240-28	28	9.3	6	20"(51cm)	10	1245-32	32	10.7	6	20"(51cm)	10
1240-32	32	10.7	6	20"(51cm)	10	1245-36	36	12.0	6	20"(51cm)	10
1240-36	36	12.0	6	20"(51cm)	10	1245-40	40	13.3	6	20"(51cm)	10
1240-40	40	13.3	6	20"(51cm)	10						

102 E. KEEFE AVENUE • MILWAUKEE, WISCONSIN 53212 U.S.A.

Tel 414/265•7620 • Fax 414/265•7628

e-mail: medo@medovations.com

1•800•558•6408