



Bougie Trade-in Policy and Recycling Service



Save 10% When You Purchase New Bougies.


Is it time to replace your esophageal bougies or wire-guided dilators?

Replace bougies when you notice:

- End cap separating from bougie
- Cuts, cracks or nicks
- Flaking or peeling of the outer coating
- Faded or missing markings
- Expired or missing expiration dates

Replace wire-guided dilators when you notice:

- Loss of flexibility
- Cuts, cracks or nicks
- Stretching in the inner lumen
- Faded or missing markings
- Curved dilator or curled tapered tip
- Expired



Trade in Hurst or Maloney tip bougies—from any manufacturer—and save 10% when you purchase a replacement from Medovations.

Trade-in Policy

- 10% discount on the purchase of new bougies
- For all Maloney and Hurst Bougies from any manufacturer
- Valid on a unit-for-unit basis

Recycling Service

- NO disposal costs
- NO worries
- NO hassles

How Do I Place My Order?

1. Contact Customer Service at 800-558-6408 or via email at orders@medovations.com to order your new bougies.
2. Communicate that your facility is replacing expired bougies and Medovations will invoice you at 10% below the selling price for each bougie being replaced.
3. If your facility is replacing mercury-filled bougies (orange-colored), inform your Product Specialist that you will need recycling kits. Mercury-free (tungsten-filled) bougies should be thrown away at your facility.
4. Recycling instructions for your expired mercury-filled bougies will be provided to you when your purchase order is placed.

Learn more from your personal product representative at **800-558-6408** or **Medovations.com**



Medovations, Inc. | 102 East Keefe Avenue | Milwaukee, Wisconsin 53212 USA
800.558.6408 or 414.265.7620 | Fax 414.265.7628
Orders: orders@medovations.com | General Inquiries: medo@medovations.com