

## Product Specifications

### Cables for use with Circon-ACMI/ Olympus and other light sources

#### 40 FR

Part No. 1450-40, InnerVision™ Fiber Optic Cable, 40 FR, Circon-ACMI/Olympus Connector, Clean, non-sterile.

#### 50 FR

Part No. 1450-50, InnerVision™ Fiber Optic Cable, 50 FR, Circon-ACMI/Olympus Connector, Clean, non-sterile.

#### 56 FR

Part No. 1450-56, InnerVision™ Fiber Optic Cable, 56 FR, Circon-ACMI/Olympus Connector, Clean, non-sterile.

### Wolf/Dyonics Adapter for Circon-ACMI/Olympus light source

Part No. 1450-02, Wolf/Dyonics Adapter for InnerVision™ Cables with a Circon-ACMI/Olympus Connector. Clean, non-sterile.

### Cables for use with Storz light source

#### 40 FR

Part No. 1455-40, InnerVision™ Fiber Optic Cable, 40 FR, Storz Connector, Clean, non-sterile.

#### 50 FR

Part No. 1455-50, InnerVision™ Fiber Optic Cable, 50 FR, Storz Connector, Clean, non-sterile.

#### 56 FR

Part No. 1455-56, InnerVision™ Fiber Optic Cable, 56 FR, Storz Connector, Clean, non-sterile.

### InnerVision™ Detachable Tips

#### 40 FR

Part No. 1465-40, InnerVision™ Detachable Tips, 40 FR, Blunt tip, box of 5, sterile.

#### 50 FR

Part No. 1465-50, InnerVision™ Detachable Tips, 50 FR, Tapered tip, box of 5, sterile.

#### 56 FR

Part No. 1465-56, InnerVision™ Detachable Tips, 56 FR, Tapered tip, box of 5, sterile.

### Light Source Compatibility Guide

InnerVision™ Fiber Optic Cable with Circon-ACMI/Olympus Connector is compatible with:

ACMI	LUXTEC
AFI	OLYMPUS
ASPEN	PILLING
CIRCON	STRYKER
CUDA	V. MUELLER
EDER	WEHMER
KLI	ZIMMER

InnerVision™ Fiber Optic Cable with Circon-ACMI/Olympus Connector and Wolf/Dyonics Adapter is compatible with:

WOLF  
DYONICS

InnerVision™ Fiber Optic Cable with Storz Connector is compatible with:

KARL STORZ (5.0 Maximum Bundle Diameter)

*Product is shipped FOB Milwaukee, WI U.S.A. Freight charges are prepaid and added to the invoice. All invoices are due upon receipt.*

### InnerVision™ Lighted Bougie

medovations, inc.  
102 east keefe avenue  
milwaukee, wi 53212 usa  
800.558.6408  
414.265.7620  
fax 414.265.7628  
medo@medovations.com  
www.medovations.com

FMSM-026 Issue

  
medovations