

CleanFreak® Endoscope Cleaning Sponges

The pre-cleaning kit holds an enzyme-imbibed sponge in a water basin.



Introducing the CleanFreak® Endoscope Cleaning Sponge Line.

At Medovations, we look for ways to make your work easier. So we've added smart Cleaning Sponges to our line in convenient designs that ensure quick and confident pre-cleaning.

Dry Imbibed Sponges. Just add water.


Infused with a water-activated enzymatic cleaner, there's no need to inventory additional solutions. There's no mess with Medovations sponges compared to wet sponge packs.

Multi-enzyme detergent.


To remove bio-burden, our imbibed sponges contain protease, amylase and lipase as compared to other dual-enzyme competitors. Our dry imbibed sponges have a low-foam formula for endoscope reprocessing.

Durable cleaning solutions that fit the way you work.

If you prefer to use a specific cleaner, choose non-imbibed sponges in cylindrical and flat textured designs. Whichever way you prefer to work, our economical sponges will help you get the job done.



The contoured design fully surrounds the endoscope for thorough cleaning.



Textured sponges conform to any product shape.

CleanFreak® Endoscope Cleaning Sponges

Part No.	Description	Qty/Case	Color
1201-01	Endoscope Sponge—Textured	200	Blue
1201-02	Endoscope Sponge—Cylindrical	50	Yellow
1201-03	Dry Imbibed Endoscope Sponge—Cylindrical	50	Yellow
1201-20	Pre-Cleaning Kit—Imbibed sponge, water basin	24	Sponge—Yellow Basin—Purple

Single-use, non-sterile. Made in the USA.

Learn more from your personal product representative at **800-558-6408** or **Medovations.com**



Medovations, Inc. | 102 East Keefe Avenue | Milwaukee, Wisconsin 53212 USA
 800.558.6408 or 414.265.7620 | Fax 414.265.7628
 Orders: orders@medovations.com | General Inquiries: medo@medovations.com